

Sterile Urethral Catheters

NANCHANG YILI MEDICAL INSTRUMENT CO., LTD

Office: 101# Building, Greenland International Exp Center, Jiulong Road 1177, Honggutan District, Nanchang, Jiangxi, China

Factory: Jinxian Medical Industrial Park, Nanchang, Jiangxi, China

Tel: +86-791-88252326 Fax: +86-791-88251326 Post Code: 330046


E-mail: sales@yilimed.com

www.yili-medical.com



HISTORY OF YILI MEDICAL

- ----- 2010 Established Medical Research and Development Institution
- ----- 2011 Established YiLi Medical Sales Department in Hong Kong
- ----- 2012 Established Nanchang YiLi Medical Instrument Co,Ltd in Nanchang City
- ----- 2013 Established first factory in Jinxian town of Nanchang
- ----- 2014 Got ISO and CE certificated, the same year dispatch a sales team to Mexico and Turkey for marketing survey
- ----- 2016 Invested other two family own factories to expand products range and also registered EasyThru Brand
- ----- 2017 Signed distributors in Turkey, Korea, Germany
- ----- 2019 Obtained tenders in Qatar, UK
- ----- 2020 Registered disposable face mask in Germany
- ----- 2021 Entered online exhibition due to COVID-2019
- ----- 2023 Registered in U.S. FDA




■ Sterile Urethral Catheters

- 1.Usage: Easythru Single Use Latex Urethral Catheter is used to empty the bladder and collect urine in a drainage bag.
- 2.Material: base material is 100% high quality latex imported from Malaysia, coated material is medical grade silicon to reducing urethral irritation and allergic reaction.
- 3.Rubber valve and plastic valve for choose. The rubber valve has our unique design to ensuring leakage rate 40-60% less than others in the market. Plastic valve design allows the caregiver to operate by one hand.
- 4.Stable flow in both drainage lumen and balloon inflation lumen
- 5.Symmetry balloon after inflating and no burst balloon according to test.
- 6.Smooth and symmetry murphy eyelets.
- 7.Certification: ISO, CE, Free Sale Certificate, cytotoxicity test report.
- 8.Hydrophilic coating is available.




EasyThru®



STRAIGHT TIP




TIEMANN TIP



■ Sterile Urethral Catheters

-Latex foley catheter siliconized

It is an indwelling catheter fasten placement by a balloon and used to empty the bladder. The third way is used for irrigation or medicine injection.



(Pediatric with stylet)



(Adult)



3 Way Rubber Valve

2Way Rubber Valve

3 Way Plastic Valve

2Way Plastic Valve


Type		Code		Size (FR)	Length (MM)	Balloon (ML)	Color
		Rubber Valve	Plastic Valve				
2 way	Pediatric with stylet	YL-LFC2R-06	YL-LFC2-06	6	270-280	3, 5, 3-5	Pink
		YL-LFC2R-08	YL-LFC2-08	8	270-280	3, 5, 3-5	Black
		YL-LFC2R-10	YL-LFC2-10	10	270-280	3, 5, 3-5	Grey
	Adult	YL-LFC2R-12	YL-LFC2-12	12	390-400	5-10, 5-15, 30, 30-50	White
		YL-LFC2R-14	YL-LFC2-14	14	390-400	5-10, 5-15, 30, 30-50	Green
		YL-LFC2R-16	YL-LFC2-16	16	390-400	5-10, 5-15, 30, 30-50	Orange
		YL-LFC2R-18	YL-LFC2-18	18	390-400	5-10, 5-15, 30, 30-50	Red
		YL-LFC2R-20	YL-LFC2-20	20	390-400	5-10, 5-15, 30, 30-50	Yellow
		YL-LFC2R-22	YL-LFC2-22	22	390-400	5-10, 5-15, 30, 30-50	Purple
		YL-LFC2R-24	YL-LFC2-24	24	390-400	5-10, 5-15, 30, 30-50	Blue
		YL-LFC2R-26	YL-LFC2-26	26	390-400	30, 30-50	Pink
		YL-LFC2R-28	YL-LFC2-28	28	390-400	30, 30-50	Brown
		YL-LFC2R-30	YL-LFC2-30	30	390-400	30, 30-50	Grey
		3 way	Adult	YL-LFC3R-16	YL-LFC3-16	16	390-400
YL-LFC3R-18	YL-LFC3-18			18	390-400	30, 30-50	Red
YL-LFC3R-20	YL-LFC3-20			20	390-400	30, 30-50	Yellow
YL-LFC3R-22	YL-LFC3-22			22	390-400	30, 30-50	Purple
YL-LFC3R-24	YL-LFC3-24			24	390-400	30, 30-50	Blue
YL-LFC3R-26	YL-LFC3-26			26	390-400	30, 30-50	Pink
YL-LFC3R-28	YL-LFC3-28			28	390-400	30, 30-50	Brown

■ Sterile Urethral Catheters


- Tiemann Tip Latex Foley Catheter Silconized

EasyThru®

Usage: It is an indwelling catheter fasten placement by a balloon and used to empty the bladder. If you have an enlarged prostate gland, stricture or scarring, the tiemann tip help to overcome obstructive difficulties when catheterization.



(TIEMANN TIP)



Type	Code		Size (FR)	Length (MM)	Balloon (ML)	Color	
	Rubber Valve	Plastic Valve					
TIEMANN TIP	2 way	YL-TTLFC2R-12	YL-TTLFC2-12	12	390-400	5-10, 5-15, 30, 30-50	White
		YL-TTLFC2R-14	YL-TTLFC2-14	14	390-400	5-10, 5-15, 30, 30-50	Green
		YL-TTLFC2R-16	YL-TTLFC2-16	16	390-400	5-10, 5-15, 30, 30-50	Orange
		YL-TTLFC2R-18	YL-TTLFC2-18	18	390-400	5-10, 5-15, 30, 30-50	Red
		YL-TTLFC2R-20	YL-TTLFC2-20	20	390-400	5-10, 5-15, 30, 30-50	Yellow
		YL-TTLFC2R-22	YL-TTLFC2-22	22	390-400	5-10, 5-15, 30, 30-50	Purple
		YL-TTLFC2R-24	YL-TTLFC2-24	24	390-400	5-10, 5-15, 30, 30-50	Blue
	3 way	YL-TTLFC3R-16	YL-TTLFC3-16	16	390-400	5-10, 5-15, 30, 30-50	Orange
		YL-TTLFC3R-18	YL-TTLFC3-18	18	390-400	5-10, 5-15, 30, 30-50	Red
		YL-TTLFC3R-20	YL-TTLFC3-20	20	390-400	5-10, 5-15, 30, 30-50	Yellow
		YL-TTLFC3R-22	YL-TTLFC3-22	22	390-400	5-10, 5-15, 30, 30-50	Purple
		YL-TTLFC3R-24	YL-TTLFC3-24	24	390-400	5-10, 5-15, 30, 30-50	Blue
		YL-TTLFC3R-26	YL-TTLFC3-26	26	390-400	5-10, 5-15, 30, 30-50	Pink
		YL-TTLFC3R-28	YL-TTLFC3-28	28	390-400	5-10, 5-15, 30, 30-50	Brown


■ Sterile Urethral Catheters

- Latex Hematuria Catheter

1. Usage: after prostatectomy, prostate cancer surgery, bladder cancer resection and other surgery operation, the catheter is used for quick compression hemostasis by big balloon, quick suction away blood clot, drainage hard substance, bladder irrigation, urine drainage, drug injection.
2. Super flow rate, effectively prevent blood clotting and suction away blood clot and tissue after operation.
3. Balloon injection volume is up to 80ml, blasting limit is up to 250ml
4. Dufour tip and couvelaire tip for choose.
5. Reinforced tube wall for choose. It could anti-collapse during clot suction under negative pressure 80 mpa.
6. With exclusive patent.



Hematuria Catheter balloon size:50-80ml



Standard balloon size: 30ml

EasyThru®

Type	Code		Size (FR)	Length (MM)	Balloon (ML)	Color
	Normal Tube Wall	Reinforced Tube Wall				
Dufour Tip	YL-LFCHD3-18	YL-LFCHDR3-18	18	390-400	50-80	Red
	YL-LFCHD3-20	YL-LFCHDR3-20	20	390-400	50-80	Yellow
	YL-LFCHD3-22	YL-LFCHDR3-22	22	390-400	50-80	Purple
	YL-LFCHD3-24	YL-LFCHDR3-24	24	390-400	50-80	Blue
Couvelaire Tip	YL-LFCHC3-18	YL-LFCHC3R-18	18	390-400	50-80	Red
	YL-LFCHC3-20	YL-LFCHC3R-20	20	390-400	50-80	Yellow
	YL-LFCHC3-22	YL-LFCHC3R-22	22	390-400	50-80	Purple
	YL-LFCHC3-24	YL-LFCHC3R-24	24	390-400	50-80	Blue

■ Sterile Urethral Catheters


- Double Balloon 3 Way Latex Urethral Catheter

Usage: the front 30 ml balloon is to fasten placement and the back 50 ml balloon is to stop the blood by compression.



EasyThru®



A(distance between two balloon 0.5 cm)



B(distance between two balloon 3.5 cm)




Type	Code	Size (FR)	Length (MM)	Balloon (ML)	Color
A (distance between two balloon 0.5 cm)	YL-BD3A-18	18	390-400	30 50	Red
	YL-BD3A-20	20	390-400	30 50	Yellow
	YL-BD3A-22	22	390-400	30 50	Purple
	YL-BD3A-24	24	390-400	30 50	Blue
	YL-BD3A-26	26	390-400	30 50	Pink
B (distance between two balloon 3.5 cm)	YL-BD3B-18	18	390-400	30 50	Red
	YL-BD3B-20	20	390-400	30 50	Yellow
	YL-BD3B-22	22	390-400	30 50	Purple
	YL-BD3B-24	24	390-400	30 50	Blue
	YL-BD3B-26	26	390-400	30 50	Pink


■ Sterile Urethral Catheters

- LATEX MALECOT CATHETER

1. Usage: Latex Malecot Catheter, also called Latex Pezzer Catheter, is used for drainage of urine directly from bladder or kidney.
2. Mushroom tip with three holes.



EasyThru[®]



Code	Size (FR)	O.D (MM)	Length (MM)
YL-LMC-12	12	4.0	400
YL-LMC-14	14	4.7	400
YL-LMC-16	16	5.3	400
YL-LMC-18	18	6.0	400
YL-LMC-20	20	6.7	400
YL-LMC-22	22	7.3	400
YL-LMC-24	24	8.0	400
YL-LMC-26	26	8.7	400
YL-LMC-28	28	9.3	400
YL-LMC-30	30	10.0	400
YL-LMC-32	32	10.7	400
YL-LMC-34	34	11.3	400
YL-LMC-36	36	12.0	400

■ Sterile Urethral Catheters

- Packaging


EasyThru®



(With inner sheath)



(Blister individual)



(Box)

**Packing: Blister individual ,10pcs/box,500pcs/ctn.
CTN : 52*41*45cm.
135000PCS/20GP, 291000PCS/40GP.**

